

1 IMPORTS MENU

The **IMPORTS MENU** offers options to import (upload) files directly to the state registry. The options that display depend on your user access and permissions. The following options are available:

- **DTT/HL7** -- This option allows users to send patient information files that originate in the facility's patient management system (or other source that matches the import type requirements). This includes Remote Clinic contraindications, and ENCD Data.
 - eNCD is used to notify the state when a person under the age of 18 has died. We can accept this data so that we can stop mailings from being sent to deceased children. The ENCD file can be loaded through the DTT and run through the normal de-duplication logic. All records coming from this file will be marked with the "inactive status" of deceased. If a record already exists for the patient, that record will be updated. If a record does NOT exist, a new record will be created.
- **VacMan Upload** This option allows users to upload patient information from the Vaccine Manager, referred to as VacMan.
- **HL7 Realtime** This option allows users to upload patient information via HL7 in Realtime.

There are different file specifications depending on the billing system.

Once the file is imported by the facility, the state will import it into the registry.



Using a scheduler such as the Microsoft Scheduler can also setup automatic running/processing of files or reports such as the DTT Imports.

When importing batch vaccinations, the vaccination that ends up in the master record is the one that contains the most information; i.e., specific vaccine type, administered shot, lot number, combination vaccine. A procedure is used to "score" the immunization to determine if the incoming or existing shot is the "winner" for the master record. The score is determined as follows:

- Administered shot (historical is null) 1
- Specific Vaccine type 2
- Has Lot Number 3
- Combination Vaccine 1
- Total Possible Score

ABOUT DTT/HL7 IMPORTS

The following can be performed from the "DTT/HL7 Import" window:

- Import (an existing Data Profile)
- Define (a New) Data Profile
- Define a new Profile with Pre-Defined Fields
- Edit an Existing Data Profile
- Print a Data Profile Report
- View Import Logs
- The maximum Import file size is 100 MB.

The **DTT** import process will correctly identify the records with the following common problems from the ACS file.



- 0Y (zero + Y) can become OV (letters O + V).
 - Any letter O in the Mail Code of the Keyline should be treated as a Zero
 - Any letter V in the Mail Code of the Keyline should be treated as a letter Y.

Language Code can be incorrect.

- If the language code contains an S assume it is Spanish. Otherwise, assume it is English.
- If a mail code starts with I assume it is INPKT.

As of version 3.0, a state configurable option can be set to prevents short names from overwriting full names during a batch DTT/HL7 import. For example, if Chris were imported, the short name of Chris would not overwrite the full name of Christopher.

Continue to the section corresponding to the action to perform.

IMPORTING DTT OR HL7 FILES

This option assumes a data profile has already been defined and you are ready to upload the file to the State.


Note:

In order to make file uploads more efficient, DTT import files may be zipped. Multiple files may be included in a single ZIP file, as long as all files match the profile being used for the import.

To import a DTT or HL7 file, perform the following:


1. From the **IMPORTS** menu point and click on the **DTT/HL7** option. The "DTT/HL7 Import" window appears.





- 2. In the **SELECT IMPORT PROFILE** field, click the drop-down arrow to view a list of valid import profiles and click one to select it. A few possible examples are:
 - PCIC Demographic
 - PCIC Facility
 - PCIC Immunization
 - PCIC Physician
 - Remote Clinic
 - Gateway Message Format
 - a) Can contain up to four unique files (Facility (F), Physician (Y), Demographic (D), and Immunizations (I)).
- **3.** For states using Data ID Mapping, the user will need to select the Import Profile AND a Facility/Physician Mapping for the IRMS from a drop-down menu.
- 4. In the **SELECT FILE** field, click the **BROWSE** button and select the file to import. The "Choose File" window appears.







- **5.** Point and double-click the file name to select it or click the **OPEN** button. The "DTT/HL7 Import" window reappears with the "chosen" file name.
- 6. In the "CLICK HERE TO PREVENT REFRESHING PATIENT DEMOGRAPHIC DATA" checkbox, click it to:
 - Insert patient demographic records as "historical" since they already exist, the demographic information in the master patient record will NOT get updated.
- 7. If the patient doesn't exist in the registry, they're added as usual
- 8. In the "CLICK HERE TO DECREMENT VACCINE INVENTORY" checkbox, click it to decrease the provider's vaccine inventory; otherwise, vaccinations are not decremented in the provider's vaccine inventory.
 - If a vaccination lot number is sent through the DTT and matches a lot number in the provider's inventory, decrement the number of vaccinations in the provider's inventory.
 - If a vaccination lot number is sent through the DTT and **DOES NOT MATCH** a lot number in the provider's inventory, do not decrement. The



- vaccination will be rejected and noted on the DTT Error Report.
- If a vaccination is sent through the DTT and the vaccination lot number is inactive, the vaccination will be rejected and noted on the DTT Error Report.
- **9.** Point and click on one of the available buttons:
 - **IMPORT** button to send the file to be processed by the State Registry staff.
 - This may take some time to process depending on the size of the file. The browser may time out when processing lengthy imports. Periodically, check the "Import Log" to ensure the import processed properly.
 - **CLEAR** button to clear the fields in order to reenter the criteria.

CREATE A NEW DATA PROFILE OR EDIT EXISTING PROFILE

This button is used when you want to create a new data profile, create a profile with "pre-defined fields," or edit an existing profile.

1. From the "DTT/HL7 Import" window, click the **DEFINE DATA PROFILE** button. The "DTT Profile Creator –
Add/Select" window appears.

Note: Figure 1-3: DTT Profile Creator – Add/Select Profile




If you are creating a NEW PROFILE, continue to step two; otherwise if editing an EXISTING PROFILE, click the drop-down arrow in the EXISTING PROFILES field and select a profile name.

6



- 2. If you are creating a **NEW DATA PROFILE**, type the new profile name in the **NEW PROFILE** field. If editing an **EXISTING PROFILE**, click the drop-down arrow for a list of existing profiles.
- **3.** Click one of the available buttons:
 - **CANCEL** to not create/edit the profile and exit to the "DTT/HL7 Import" window.
 - **CLEAR** to erase the fields and re-enter profile names/choices.
 - **CONTINUE** to access the file format selection. The "DTT Profile Creator Select Format" window appears.

Note: Figure 1-4: DTT Profile Creator – Select Format



- 4. Type an example of the incoming data in the **SAMPLE** LINE field.
 - For demonstration purposes (shown below), we have a delimited file using the pipe symbol for our delimiter. The sample line will contain a medical record number, first name, last name, date of birth, and gender.
- **5.** Indicate the incoming file's field separator type by clicking one of the radio buttons:
 - **TAB Delimited** indicates the data fields are separated using TABS.
 - **DELIMITED** indicates a special symbol is used to separate the data fields.
 - a) Type the symbol used in the **DELIMITER** field. For example, type the pipe symbol (|).



- When using Full Names (Patient, Guardian, Facility Contact, or Physician), you CANNOT use a comma as the delimiter symbol.
- **FIXED WIDTH** indicates the data fields are fixed in width.
 - a) Type a number for each field's length, separating each with a comma.

Note: Figure 1-5: DTT Profile Creator – Select Format with Sample Line



6. The fields and their descriptions are listed in the table:

Note: Table 1-1: DTT Profile Creator – Select Format Field Descriptions

COLUMN	DESCRIPTION		
SAMPLE LINE	A sample row of the data format being imported.		
	Note: It is recommended to copy a sample line from the actual file being imported, and then paste it into the Sample Line to alleviate errors.		
TAB DELIMITED	Indicates the data file being imported uses the TAB symbol/character as the separator between the fields.		
	You MUST copy and paste from the file being imported.		
	Be sure to click the Tab Delimited radio button.		
DELIMITED	• Indicates the data file being imported uses a special symbol between fields. For example, the pipe symbol (). Copy/paste, or type the sample line and indicate what the delimiter character is in the DELIMITER field.		



COLUMN	DESCRIPTION
	• Note: When using FULL NAMES (Patient Full Name, Guardian Full Name, Facility Contact Full Name, or Physician Full Name), the COMMA symbol CANNOT be used as the separator character since it is used to define the Data Format.
	Be sure to click the Delimited radio button.
FIXED WIDTH	• Indicates the data file being imported uses a fixed width between fields. For example, if the first field uses 25 characters, and the second field uses 10, your entry in the Field Length would be 25,10
	Be sure to click the Fixed Width radio button.

7. Click one of the available buttons:

- **CANCEL** to return to the "DTT Profile Creator Add/Select Profile" window.
- **CLEAR** to erase the data, remain on the window to re-enter new data.
- **EVALUATE SAMPLE LINE** to continue and access the "DTT Profile Creator—Assign Fields" window.



Note: Figure 1-6: DTT Profile Creator – Assign Fields

DTT Profile Creator - Assign Fields						
Sample Field	Field Definition		Requires Additional Field Defs		Data Format	
12345	Medical Record Number	¥		-select	t-	~
Sherri	Patient First Name	~		-select	t-	~
Bull	Patient Last Name	~		-select	t-	~
06031960	Patient DOB	~		MM/DE	DWW	~
F	Patient Mailing Address1	٧		-select	t-	~
		Ca	ncel R	eset	Submit Definitions	

8. The columns and their descriptions are listed in the table:

Note: Table 1-2: DTT Profile Creator – Assign Fields Descriptions

	Difficulty Presign Fields Descriptions
COLUMN	DESCRIPTION
SAMPLE FIELD	A sample row of data that is incoming. This is the sample line that was entered on the "DTT Profile Creator – Select Format" window and cannot be changed here.
FIELD DEFINITION	This column is used to indicate IWEB's description of the Sample (Line) Field so it can be identified. Click on the drop-down menu arrow to view/select a valid choice. Examples are:
	Patient Full Name
	Guardian Full Name
	Facility Contact Full Name
	Physician Full Name
	Note: It is recommended to choose all of the Field Definitions before choosing other information.
REQUIRES ADDITIONAL FIELD DEFS	This column is used to indicate that a field will also be used as an additional field in IWEB. For example, if a medical record number is entered, this could also be mapped as a Patient Birth File Number.
	Another example would be for multiple race codes (White and African American) would need to be mapped to indicate there are multiple codes.
	Be sure to click the checkbox if additional mapping needs to be defined.




COLUMN	DESCRIPTION
DATA FORMAT	This column is used to indicate the format of the incoming field. Click on the dropdown menu arrow to view/select the appropriate format. Examples are:
	Dates
	First Middle Last Suffix
	First Middle Last, Suffix
	Last Suffix, First Middle
	Last, First, Middle, Suffix
	Last, First Middle Suffix
	Both Middle (name) and Suffix are optional. If the format selected does not have commas to differentiate the location of each, and only one is passed, the value will be considered a suffix if it is one of the following:
	• II, III, IV, SR, JR (this list includes the spelling variations, Junior, Senior, 2 nd , etc.)
	Validation has been added to name formats. If a name format is selected and the data passed is not in the correct format (i.e., name is supposed to have commas and it does not), the user will receive an error message viewable via the DTT Import Log.

- 9. Using the sample field, choose the appropriate **FIELD DEFINITION** by clicking the drop-down menu arrow to view/select a valid definition. For example, the first line in the "Sample Field" (shown above) contains the numbers 12345.
 - Click on the drop-down arrow in the "Field Definition" column, and choose "medical record number."
 - Continue defining the fields until all have been defined.

The window should now resemble the one shown below.



Figure 1-7: DTT Profile Creator - Assign Fields with Field Definitions



- **10.** Using the Field Definition column, click the drop-down menu arrow to view/select a valid **DATA FORMAT** specifically for:
 - Dates
 - First Middle Last Suffix
 - First Middle Last, Suffix
 - Last Suffix, First Middle
 - Last, First, Middle, Suffix
 - Last, First Middle Suffix
- 11. To indicate whether a field needs additional mapping (e.g., multiple race codes), click the **REQUIRES ADDITIONAL FIELD DEFS** checkbox.
- **12.** Click one of the available buttons:
 - **CANCEL** to return to the "DTT/HL7 Import" window.
 - **CLEAR** to erase and re-enter new field/date definitions.
 - **SUBMIT DEFINITIONS** to save the profile and display one of the following windows:
 - a) "DTT Profile Creator Define Multiple Value." Continue to the section titled, "Defining Multiple Values."
 - b) "DTT Profile Creator Required Mappings." Continue to the section titled, "Mapping Fields."




value." Continue to the section titled, Defining Default Values."

DEFINING MULTIPLE VALUES

The "DTT Profile Creator – Define Multiple Value" window appears when you want to map a field to have an additional value. For example, medical record number could be the same as Patient Birth File number.


Note: Figure 1-8: DTT Profile Creator – Define Multiple Value



1. From the "DTT Profile Creator – Define Multiple Value" window, click the **SELECT ARROW** button. The "DTT Profile Creator – Define Multiple Values – (Edit)" window appears:



Note: Figure 1-9: DTT Profile Creator – Define Multiple Values (Edit)



- **2.** Click the drop-down arrow to view/select an additional field for the Sample Field value.
- **3.** Click one of the available buttons:
 - **RETURN TO MENU** returns to the "DTT Profile Creator Define Multiple Value" window. If you have not clicked the Add or Remove button, no changes will be saved.
 - ADD DEFINITION saves the selection, adds the DELETE checkbox, and remains on the "DTT Profile Creator Define Multiple Values (Edit)" window.

Note: Figure 1-10: DTT Profile Creator – Define Multiple Values (Edit) with Delete Checkbox



- **REMOVE DEFINITION** deletes the additional field value after clicking in the **DELETE** checkbox and remains on the "DTT Profile Creator Define Multiple Value (Edit)" window.
- **4.** Click the **RETURN TO MENU** button. The "DTT Profile Creator Define Multiple Value" window appears.
- **5.** Click the **CONTINUE** button. The "DTT Profile Creator Required Mappings" window appears. Continue to the next section titled, "Mapping Fields" for instructions.



MAPPING FIELDS

The word "mapping" means the translation of the import file codes coinciding with the registry acceptance codes).


The "DTT Profile Creator – Required Mappings" window appears after clicking the CONTINUE button on the "DTT Profile Creator – Define Multiple Value" window.

Examples of when you might see fields listed in this window are when your profile includes the fields:

- Patient Gender Code
- Patient Ethnicity Code
- Vaccine Code
- Manufacturer Code
- County Code, etc.

If you don't see the "DTT Profile Creator – Required Mappings" window, continue to the next section.


Note: Figure 1-11: DTT Profile Creator – Required Mappings




1. Click the **HYPERLINK** (words underscored in blue). The "DTT Profile Creator – Define (name of field)" window appears.



Note: Figure 1-12: DTT Profile Creator – Define (Name of Field)



- 2. In the **CODE** column, type the corresponding code that matches the **DESCRIPTION** column.
- **3.** Click one of the available buttons:
 - **RESET** to erase any newly typed text and remain on the same window.
 - **SUBMIT DATA MAPPINGS** to save and return to the "DTT Profile Creator Required Mappings" window.




- **4.** Click one of the available buttons:
 - **CANCEL** to not save and exit to?
 - **NEXT** to save and display the "DTT Profile Creator Define Default Value" window. Continue to the next section titled, "Defining Default Values" for instructions.

DEFINING DEFAULT VALUES

The "DTT Profile Creator – Define Default Value" window appears after clicking NEXT on the "DTT Profile CREATOR – Required Mappings" window. It is used to set up defaults for specific fields. An example of use is the Facility's Address State Code.


Note: Figure 1-13: Profile Creator – Define Default Value





- 1. In the **FIELD DEFINITION** column, click the drop-down menu arrow to view/select a **FIELD**.
- 2. Press the **TAB** key to position the mouse cursor in the **DEFAULT VALUE** field.
- **3.** Type the appropriate **VALUE** in the field and then click one of the available buttons:
 - **ADD NEW DEFAULT** saves the entries and places the new default at the bottom of the window.

Figure 1-14: DTT Profile Creator – Define Default Value – Default Added




- **CLEAR** erases any selections/text entered and remains on the window.
- **REMOVE DEFAULT** Removes the default value and remains on the window.
 - This button will not appear if there aren't any defaults set.
- **SAVE PROFILE** Saves the profile and returns to the "DTT/HL7 Import" window. If you have added a new default, verify it displays in the bottom



portion of the window before saving the profile, or your defaults will not be saved.

Note: Figure 1-15: DTT/HL7 Import – Profile Saved Successfully



Note:

PRINT DATA PROFILE REPORT

This option will print an existing data profile report.

- 1. From the "DTT/HL7 Import" window, in the **SELECT IMPORT PROFILE** field, point and click the drop-down arrow to view/select a valid profile.
- 2. Click the **PRINT DATA PROFILE REPORT** button to display the report in a new browser window providing the option to print.



Figure 1-16: Data Profile Report

PCIC DEMOGRAPHIC

IRMS: 1000 **Creation Date:** Last Updated:

Format: Delimited Delimiter: |

Field order and mappings:

- · Medical Record Number
- Patient First Name
- Patient Middle Name
- Patient Last Name
- · Patient Suffix
- Patient SSN
- Patient DOB
 - o Date Format: MM/DD/YYYY
- Patient Mailing Address1
- Patient Physical Address
- Patient Address City
- · Patient Address State Code
- Patient Address Zip Code
- Patient Gender Code
- · Patient Race Code
- Mother Maiden Name
- · Guardian First Name
- Guardian Middle Name
- · Guardian Last Name
- · Guardian SSN
- Patient Phone Number Area Code
- · Patient Phone Number
- Inactive Code
- Patient Medicaid Number
- Managed Care Plan Name
- Managed Care Plan ID
- Managed Care Plan Date
 - o Date Format: MM/DD/YYYY
- Patient County FIPS CodePatient Ethnicity Code
- Patient Health District
- · Patient VFC Eligible Code
- · Facility ID
- Physician ID
- Undefined Field
- Patient Email Address
- Block Recall
- Programs List
- School District Code
- Patient Birth State
- Patient Birth Country
- Comments
- Patient Language



VIEW IMPORT LOGS (FOR ERRORS)

This option is used to view previously submitted imports and to see what the errors are so they can be corrected and then be reimported.

1. From the "DTT/HL7 Import" window, click the **VIEW IMPORT LOGS** button. The "DTT Import Logs" window appears.

Note: Figure 1-17: DTT Import Logs



2. The columns that appear are listed in the table:


Note: Table 1-3: DTT Import Logs Column Descriptions

COLUMN	DESCRIPTION	
VIEW ERRORS	Contains the SELECT ARROW to select the desired file log.	
PROFILE NAME	Name of the Data Profile used for the import.	
IMPORT FILE	Name of the Import File.	
DATE	Year, month, day, and timestamp the import log was created.	
TOTAL RECORDS	Number of records sent with the import.	
ERRORS	Number of errors within the error log. Note: If the total number of errors exceeds 100, you may want to import a smaller file by copying/pasting a subset of the larger "original" file and view its errors. When the number of errors is large, the browser can experience a "memory limitation" and you will not be able to view the errors.	
COMPLETED	This column indicates whether or not the full import file was read in during the import and will contain a Yes or No.	



3. Locate the file whose errors you want to view by clicking on the **SELECT ARROW** button (under the "View Errors" column). The "DTT Import Errors" window appears.

Note: Figure 1-18: DTT Import Errors



4. The fields/columns of information are listed in the table:

Note: Table 1-4: DTT Import Errors Field/Column Descriptions

COLUMN	DESCRIPTION				
FILENAME	• the file's	Name of the file that was imported appended with the file's type. The example reflects DTT Stuff.txt.			
LINE NUMBER	• error.	The row/line number	r in the file t	hat caused the	
DETAILS	The VIEW DETAILS hyperlink displays the "DTT Import Error Detail Log" for the specific line of the error. DIT Import Error Detail - Line2 Errors Found No valid records sent. At least one record type (Patlent, Vaccination, Facility, Physician) must have the minimum required fields. Back to Error List Back to Error List				
		Field Definition	Date Format	Field Value	
		Medical Record Number		null	
		Patient Full Name Patient DOB Patient Gender Code Patient Birth State	MMDDYYYY	AZ	
				Back to Error List	
ERRORS	• was car	Description of the er used by a blank row b		ror in this example	

- **5.** Click on an available button or hyperlink:
 - **CANCEL** to return to the "DTT Import Logs" list window.
 - **CREATE ERROR REPORT** to display a report of all the errors in a new browser window. The window provides an option to print. Close the browser window when finished viewing/printing.




Profile: TEST1
Filename: DTT Stuff bd
Import Date: 2005-05-10 09.42
Total Error Records: 2
Total Error Records: 1

Medical Record Number Patient DOB Errors

- No valid records sent. At least one record type (Patient, Vaccination, Facility, Physician) must have the minimum required fields.

• **VIEW DETAILS** – to display the "DTT Import Error Detail" window (sample shown below).

Note: Figure 1-20: DTT Import Error Detail – Line 1



6. Click the **BACK TO ERROR LIST** to return to the "DTT Import Errors List" window.

Types of Import Errors

If the Administrator has set the DTT Properties to "Error on Unspecified Vaccination Import," you will receive an error indicating, "The vaccine given for an administered vaccination cannot be unspecified."

At least one record type (Patient, Vaccination, Facility, or Physician) must have the minimum required fields or no valid records will be sent.

The field values and how an error could occur are listed in the table.



Table 1-5: DTT Import Errors

FIELD		VALUE MUST BE
FACILITY ZIP CODE	• long.	Numeric and either be 5 or 9-digits
FACILITY COUNTY CODE	•	Numeric or an error will occur.
FACILITY STATE CODE	•	2-characters long.
FACILITY HEALTH DISTRICT	•	Numeric or an error will occur.
PHYSICIAN SSN	• will oc	Numeric and be 9-digits long or an error cur.
PHYSICIAN HEALTH DISTRICT	•	Numeric or an error will occur.
CPT VACCINE CODE	•	Numeric or an error will occur.
CDC VACCINE CODE	•	Numeric or an error will occur.
SIIS VACCINE CODE	•	Numeric or an error will occur.
IMMUNIZATION DATE	•	Formatted correctly
	•	Cannot be a future date
	• date.	Cannot be before the patient's birth
HISTORICAL FLAG	•	Must have a Y or N for its value.
VACCINATION HEALTH DISTRICT	•	Numeric or an error will occur.
VIS PUBLICATION DATE1	•	Formatted correctly.
VIS PUBLICATION DATE2	•	Formatted correctly.
VIS PUBLICATION DATE3	•	Formatted correctly.
VIS PUBLICATION DATE4	•	Formatted correctly.
VIS FORM GIVEN DATE	•	Formatted correctly.
DELETION DATE	•	Formatted correctly.
CONTRAINDICATION VACCINE CODE	•	Numeric or an error will occur.
CONTRAINDICATION	•	Numeric or an error will occur.



FIELD		VALUE MUST BE
CODE		
LEAD SAMPLE DRAW	•	Formatted correctly
	•	Cannot be a future date
	• date.	Cannot be before the patient's birth
LEAD SAMPLE	•	Formatted correctly.
ANALYZED DATE	• Draw D	Cannot be before the Lead Sample Date.
	• date.	Cannot be before the patient's birth
LEAD RESULT	•	Formatted correctly.
REPORTED DATE	• date.	Cannot be before the patient's birth
PATIENT'S BIRTH DATE	•	Formatted correctly
	•	Cannot be a future date
PATIENT SSN	•	Numeric and be 9-digits long.
PATIENT ZIP CODE	• or an e	Numeric and be either 5 or 9-digits long rror will occur.
PATIENT COUNTY CODE	•	Numeric or an error will occur.
PATIENT PHONE EXTENSION	•	Numeric or an error will occur.
PATIENT HEALTH DISTRICT	•	Numeric or an error will occur.
GUARDIAN SSN	•	Numeric and be 9-digits long.
GUARDIAN WORK PHONE	•	Numeric and be 10-digits long.
PATIENT PHONE NUMBER	•	Numeric or an error will occur.
PATIENT FAX NUMBER	•	Numeric or an error will occur.
PHYSICIAN PHONE NUMBER	•	Numeric or an error will occur.
FACILITY PHONE NUMBER	•	Numeric or an error will occur.
MEDICAL RECORD NUMBER	• Vaccin	Entered for a Patient Record and a ation Record.



FIELD	VALUE MUST BE
PATIENT FIRST NAME	Entered for a Patient Record.
	• Cannot contain the characters '%', '_', or '?'.
PATIENT LAST NAME	Entered for a Patient Record.
	Cannot contain characters '%', '_', or '?'.
PATIENT BIRTH DATE	Entered for a Patient Record.
MOTHER MAIDEN NAME OR STREET ADDRESS/ZIP CODE OR	One of these fields must be entered for a Patient Record.
STREET ADDRESS/CITY	
VACCINE CODE	Entered for a Vaccination Record.
IMMUNIZATION DATE	Entered for a Vaccination Record.
FACILITY ID	Entered for a Facility Record.
PATIENT	Valid for a Lead Record
DEMOGRAPHIC RECORD	Valid for a Contraindication Record.
PHYSICIAN ID	Entered for a Physician Record.
LOT NUMBER	Value must match lot numbers in the inventory for the current IRMS and/or facility.
LOT NUMBER	Value should not be blank if decrementing from inventory.
MANUFACTURER CODE	Value should not be blank if decrementing from inventory.
BIRTH STATE CODE	Value must be 2-characters long.
PATIENT STATE CODE	Value must be 2-characters long.
CONTRAINDICATION DELETION DATE	Formatted correctly or an error will occur.
PERMANENT CONTRAINDICATION FLAG	Values must be Y or N.

VACMAN UPLOAD

This option is used to import the vaccine file. The SIIS Administrator must set the VACMAN UPLOAD PATH in the



Administration Properties file, and the USER must have VacMan Upload Permissions.

To import a VacMan file, perform the following:

1. From the **IMPORTS** menu, click on the **VACMAN UPLOAD** option. The "VacMan Upload" window appears.



Note: Figure 1-21: VacMan Upload

- **2.** Either type the file name or click the **BROWSE** button to locate your file. Double-click the file name. The filename and path will display in the **SELECT FILE** field.
- **3.** Click one of the available buttons:
 - **UPLOAD** to begin the import process to the SIIS Webserver. This may take a few moments/minutes depending on the size of the file. During this time, your browser may seem to be stuck waiting; However, it is just waiting until the file import has been completed. Please do NOT click the browser "Stop" button unless you wish to cancel the file import. Upon completion a confirmation message appears at the top of the window, "Your file has been uploaded."
 - **CLEAR** to erase the file name from the **SELECT FILE** field.

HL7 REALTIME IMPORTS


This option allows HL7 records to be imported in real time.

1. From the **IMPORTS** menu, click the **HL7 REALTIME** option. The "HL7 Realtime Server" window appears.

Note: Figure 1-22: HL7 Realtime Server


26





- **2.** This window is divided into two sections: Address and Messages Received (Batched by Days).
- **3.** If you want to view the HL7 Realtime Server's details, click on the **ADDRESS** Link. The HL7 Realtime Server window appears.

Note: Figure 1-23: HL7 Realtime Server Details



4. If you want to view the "Messages Received," click on the **SELECT** arrow (to the left of the date).



Figure 1-24: HL7 Realtime Messages Received, Rejected, and Accepted



5. To view the **MESSAGES ACCEPTED** or **RECEIVED**, click on the **SELECT** arrow. The respective window appears.

Note: Figure 1-25: Received Message Detail



6. Click the **CANCEL** button when you are finished.